### **Global CDMO For IVD Instruments**



# Venus 9000

Fully Automated Chemiluminescence Immunoassay Analyzer

Registration Certificate No: 渝械注准20232220334



### Navigate with Light | Lead a New Voyage

high stability | full automation | deep intelligence | broad scalability



# CONTENTS

\_\_\_\_\_

\_\_\_\_\_

V e n u s 9 0 0 0

### **Instrument Specifications**

### **Functional Modules**

### **Instrument Parameters**

### **Software Operations**

\_.\_..


### **Standard Accessories**

### Venus 9000 Ultra-fast Automatic Chemiluminescence Immunoassay Analyzer, Focus on Excellence

Registration Certificate No: 渝械注准20232220334

Fully automated chemiluminescence immunoassay analyzer— 900T/H Intelligent algorithms, flexible scheduling, adaptation to all complex reaction processes, ultra-long fault-free running time, precise control of the entire detection process with  $CV \leq 3\%$ , optional tips module to prevent contamination. Modular combination of analyzer+ sample input output buffer + track, compatible with biochemical, flow cytometry, and coagulation analyzers to form the automated laboratory testing system

### Fully compatible with the AE/AP//HRP principle





### • Ultra-fast with compact size


900T/H per analyzer, customizable connection with multiple analyzers, domestically exclusive

### • High reliability, deep customization

Ultra-long fault-free running time, precise control of the entire detection process with  $CV\% \leq$  3, optional tips module to prevent contamination. Compatible with biochemical, flow cytometry, and coagulation analyzers to form the automated laboratory testing system

### Floor efficiency ratio: 495T/H/m<sup>2</sup>

Dimensions (mm) : 1270\*1015\*1450 (Instrument) ; 510\*1015\*791 (Input Output Buffer)



### • Deep intelligence, full compatibility

Intelligent algorithms, flexible scheduling, adaptability to all complex reaction processes, fully compatible with the AE/AP/HRP principles

### • Comprehensive functionality, easy operation


Comprehensive functionality, covering a wide range of application scenarios 3-axis linkage for free-floating main user screen + secondary user screen. 360° rotational barcode scanning, allowing barcodes to be placed arbitrarily. Real continous reagent replacement



### **Modular combination of analyzer + sample input output buffer + track**

## **Functional modules**





Batch sample loading tray Up to 360 sample positions, automatic continuous sample loading




Sample transport system 3 tracks, high-speed transportation, sample storage buffer zone, no need for manual intervention

## **Functional Modules**



42 reagent positions, replacement at any time, continuous loading, 24-hour refrigeration at 2-8°C



Series	Venus 9000					
Detection Principle	Open platform, supports AE, AP, or HRP systems with magnetic particle					
Throughput	900T/H					
Sample Positions	360 sample positions					
Reagent Positions	42 reagent positions, with refrigeration function at 2~8°C. Replacement without an interruption in system operation. Customization available.					
Sampling	Optional steel probes or tips module					
Walk-Away Time	> 3h					
Precision	CV≤3%					
Carryover	0.1ppm					
Sample Loading	Batch sample loading by tray. Loading and unloading samples at any time. STAT assay at anytime					
LIS	Serial port, Network port, Unidirectional, Bidirectional					
Sample Containers	Standard blood collection vessel, Eppendorf tube, Hitachi cup					
Input Output Buffer Dimensions	510*1015*791 (mm)					
Instrument Dimensions	1270*1015*1450 (mm)					

# **Operation** Software

0	Status	Notification	Monitor	Maintain	Log	About			Statistics	Loaded	
	Time To		Test	Test Progress				Test Summary			
• Overview	00	:00			2	0	0	1	0	0	
Orders	HOUR	MIN			ending Ru	inning Co	mpleted	Samples	Exceptions	Abnormal	
📓 Data	M2	Idle	M1	Idle	T1		Tests			Ð	
Inventory	000	<b>.</b>	655	<b>7</b> .			0 items				
XQC/CAL	Reagents	supplies	Reagents	supplies	Output TrayI	Input TrayI	Sam	ple Type ÷ Sequen	ce No 💠 SID	Assay	
Settings		~		~							
Q	Notification	Control	Notification	Control	Output TrayII	Input Tray II					
	0		0								
	Calibration	Monitor	Calibration	Monitor	Buffer	STAT					
	00h:00m	0/0/0	00h:00m	0/0/0							
admin											
2023-11-10 14:51:57	Unregister	red	00h:00m	2/0/0	u 🌧 🌡c 🖶	* • No	new message	900		40 40 100	

### Overview

user friendly/real-time monitoring of assay tests/real-time monitoring of reagents and consumables/smart alert/ multi-language support



### Orders

Supports ordering specimen, control, and calibiration tests. With smart test data statistics, automatic test classification and display which makes it easy to search, rerun and release tests

# **Standard Accessories**





Reagent racks





Substrate bottles

Reaction vessels



Sample racks

### Reagents

#### **Tumor Marker**

CYFRA21-1 NSE AFP CEA CA125 SCCA CA15-3 CA19-9 CA72-4 HE4 PIVKA II ProGRP CA50 CA242 LN

### Fertility

FSH E2 AMH PROG β-HCG TESTO LH PRL SHBG

### Anemia

VB12 FA Ferr

### Cardiac Marker

hs-cTnI NT-proBNP BNP CK-MB MYO LP-PLA2 H-FABP

#### Thyroid

FT4 FT3 TSH TT3 TT4 TG TGAb TPOAb

### Inflammation

### Monitoring

PCT SAA IL-6

Kidney Marker Cortisol

#### Endocrine hormone GH ACTH

Stomach Marker

PG I PG II G17

### Coagulation

D-Dimer FDP PIC TAT TM tPAI.C

### Liver Marker

PIIINP CIV HA

### **Bone Metabolism**

PTH CT 25-OH VD

### Glyco Metabolism

C-P INS