

Reagents



Tumor Marker

CYFRA21-1
NSE
AFP
CEA
CA125
SCCA
CA15-3
CA19-9
CA72-4
HE4
PIVKA II
ProGRP
CA50
CA242
LN

Fertility

FSH
E2
AMH
PROG
 β -HCG
TESTO
LH
PRL
SHBG

Anemia

VB12
FA
Ferr

Cardiac Marker

hs-cTnI
NT-proBNP
BNP
CK-MB
MYO
LP-PLA2
H-FABP

Thyroid

FT4
FT3
TSH
TT3
TT4
TG
TGAb
TPOAb

Inflammation

Monitoring

PCT
SAA
IL-6

Kidney Marker

Cortisol

Endocrine Hormone

GH
ACTH

Stomach Marker

PG I
PG II
G17

Coagulation

D-Dimer
FDP
PIC
TAT
TM
tPAI.C

Liver Marker

PIIINP
CIV
HA

Bone Metabolism

PTH
CT
25-OH VD

Glyco Metabolism

C-P
INS



Global CDMO For IVD Instruments

Reliable Manufacturer

STAT anytime

Liquid level detection

Easy operation

Clot detection

Collision detection




KEYLIGHTS (SHENZHEN) | CHONGQING KEYSMILE BIOLOGICAL
TECHNOLOGY CO., LTD | TECHNOLOGY CO., LTD

Tel: 0086-19129380583
Hotline: 400-808-1186

Website: www.keylights.cn
Email: chemi@keylights.cn

KEYLIGHTS (Shenzhen) | CHONGQING KEYSMILE
TECHNOLOGY CO., LTD | BIOLOGICAL TECHNOLOGY CO., LTD



Basic Parameters

Series	Venus 9000	Venus 500	Venus 300	Venus 100	SMART 6500	SMART 500	KCARD 100
Sample positions	360	120	56	40	96	24	30
Reagent positions	42	32	32	15	24	15	5*50
Vessels	2000+	1000	1000	200	1000	130	/
Dimensions(mm)	1270*1015*1450 (Instrument) 510*1015*791 (Input Output Buffer)	1020*750*680	1020*750*680	680*535*580	1286*720*700	506*670*678	765*470*510